


Product Datasheet

Rabbit Endothelin 1/EDN1 ELISA Kit (orb546479)

Description	Rabbit Endothelin 1/EDN1 ELISA Kit
Reactivity	Rabbit
Range	3.9 pg/ml - 250 pg/ml
Tested Applications	ELISA
Immunogen	A peptide reflecting the target protein in rabbit. Sequence: C53-W73
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum and plasma (heparin, EDTA).
MW	23135 MW
Uniprot ID	P29560
Sensitivity	0.5 pg/ml
Expiration Date	Please enquire.



Rabbit
Endothelin
PicoKine
ELISA Kit
sta...