


Product Datasheet

Bovine TGF-Beta 2 ELISA Kit (orb546481)

Description	Bovine TGF-Beta 2 ELISA Kit
Reactivity	Bovine
Range	31.2 pg/ml - 2,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: A303-S414
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum and plasma (heparin, EDTA, citrate).
MW	47748 MW
Uniprot ID	P21214
Sensitivity	10 pg/ml
Expiration Date	Please enquire.



Bovine TGF-
beta 2
PicoKine
ELISA Kit
sta...