


Product Datasheet

Horse equine TGF-Beta 3 ELISA Kit (orb546492)

Description	Horse equine TGF-Beta 3 ELISA Kit
Reactivity	Equine
Range	31.2 pg/ml - 2,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: A301-S412
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum, plasma (EDTA) and urine.
MW	47328 MW
Uniprot ID	F7DR45
Sensitivity	10 pg/ml
Expiration Date	Please enquire.



Horse
equine TGF-
beta 3
PicoKine
ELISA K...