


Product Datasheet

Human Myostatin/GDF8 ELISA Kit (orb546495)

Description	Human Myostatin/GDF8 ELISA Kit
Reactivity	Human
Range	31.2 pg/ml - 2,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: D267-S375
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, tissue, serum and plasma (heparin or EDTA)
Uniprot ID	O14793
Sensitivity	10 pg/ml
Dilution Range	Dilution ratio of 1:1, concentration in serum and plasma is around 500 pg/ml.
Expiration Date	Please enquire.



Human
Myostatin/GDF8
PicoKine ELISA
Kit ...