


Product Datasheet

Rat Myostatin/GDF8 ELISA Kit (orb546497)

Description	Rat Myostatin/GDF8 ELISA Kit
Reactivity	Rat
Range	31.2 pg/ml - 2,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: D268-S376
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, tissue, serum and plasma (heparin or EDTA)
Uniprot ID	O35312
Sensitivity	10 pg/ml
Dilution Range	Dilution ratio of 1:10, concentration in serum and plasma is around 5 ng/ml.
Expiration Date	Please enquire.



Rat
Myostatin/GDF8
PicoKine ELISA
Kit St...