


Product Datasheet

Chinese Hamster CXCL14 ELISA Kit (orb546500)

Description	Chinese Hamster CXCL14 ELISA Kit
Reactivity	Hamster
Range	62.5 pg/ml - 4,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: E.coli; Immunogen sequence: S35-E111
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum and plasma (heparin, EDTA).
MW	13078 MW
Uniprot ID	095715
Sensitivity	10 pg/ml
Expiration Date	Please enquire.



Hamster
CXCL14
PicoKine
ELISA Kit
standa...