


Product Datasheet

Chinese Hamster BMP-7 ELISA Kit (orb546504)

Description	Chinese Hamster BMP-7 ELISA Kit
Reactivity	Hamster
Range	31.2 pg/ml - 2,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: CHO; Immunogen sequence: S292-H430
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	bone tissue, cell culture supernatants, serum and plasma (heparin, EDTA).
MW	49199 MW
Uniprot ID	P23359
Sensitivity	10 pg/ml
Expiration Date	Please enquire.



Hamster
BMP-7
PicoKine
ELISA Kit
standar...