


Product Datasheet

Chinese Hamster GDF5 ELISA Kit (orb546510)

Description	Chinese Hamster GDF5 ELISA Kit
Reactivity	Hamster
Range	31.2 pg/ml - 2,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: E.coli; Immunogen sequence: A376-R495
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum and plasma (heparin, EDTA).
MW	7491 MW
Uniprot ID	Q8K4X5
Sensitivity	10 pg/ml
Expiration Date	Please enquire.



Hamster
GDF5
PicoKine
ELISA Kit
standard...