


Product Datasheet

Rabbit GREM2 ELISA Kit (orb546515)

Description	Rabbit GREM2 ELISA Kit
Reactivity	Rabbit
Range	62.5 pg/ml - 4,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: E.coli; Immunogen sequence: R22-Q168
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum and plasma (heparin, EDTA).
MW	19320 MW
Uniprot ID	Q9H772
Sensitivity	12 pg/ml
Expiration Date	Please enquire.



Rabbit
GREM2
PicoKine
ELISA Kit
standard...