


Product Datasheet

Rat Serpin A12 ELISA Kit (orb546520)

Description	Rat Serpin A12 ELISA Kit
Reactivity	Rat
Range	15.6 pg/ml - 1,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: L21-P411
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum and plasma (heparin, EDTA).
Uniprot ID	Q8R4Z1
Sensitivity	10 pg/ml
Expiration Date	Please enquire.



Rat Serpin
A12 PicoKine
ELISA Kit
Standa...