


Product Datasheet

Human IL-1RAcP/IL-1R3/IL1RAP ELISA Kit (orb546524)

Description	Human IL-1RAcP/IL-1R3/IL1RAP ELISA Kit
Reactivity	Human
Range	312 pg/ml - 20,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: S21-T367
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum and plasma (heparin, EDTA).
Uniprot ID	Q9NPH3
Sensitivity	50 pg/ml
Dilution Range	Dilution ratio of 1:50, concentration in serum and plasma is around 200 ng/ml.
Expiration Date	Please enquire.



Human IL-1RAcP/IL-1R3/IL1RAP PicoKine EL...