


Product Datasheet

Human Omentin/intelectin-1 ELISA Kit (orb546539)

Description	Human Omentin/intelectin-1 ELISA Kit
Reactivity	Human
Range	62.5 pg/ml - 4,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: T19-S298
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum and plasma (heparin, EDTA, citrate).
Uniprot ID	Q8WWA0
Sensitivity	10 pg/ml
Dilution Range	For serum, dilution ratio is 1:1 and concentration is around 100 pg/ml. For plasma, dilution ratio is 1:5 and concentration is around 2.5 ng/ml.
Expiration Date	Please enquire.



Human
Omentin/intelectin-
1 PicoKine ELIS...