


Product Datasheet

Mouse SEZ6L ELISA Kit (orb546551)

Description	Mouse SEZ6L ELISA Kit
Reactivity	Mouse
Range	78 pg/ml - 5,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: L29-A958
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum and plasma (heparin).
Uniprot ID	Q6P1D5
Sensitivity	50 pg/ml
Dilution Range	Dilution ratio of 1:1, concentration in serum and plasma is 2-3 ng/ml.
Expiration Date	Please enquire.



Mouse
SEZ6L
PicoKine
ELISA Kit
Standard ...