


Product Datasheet

Human SLITRK6 ELISA Kit (orb546556)

Description	Human SLITRK6 ELISA Kit
Reactivity	Human
Range	15.6 pg/ml - 1,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: S27-S608
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum and plasma (heparin).
Uniprot ID	Q9H5Y7
Sensitivity	10 pg/ml
Dilution Range	Dilution ratio of 1:1, concentration in serum and plasma is less than the lowest standard, 15.6 pg/ml.
Expiration Date	Please enquire.



Human
SLITRK6
PicoKine
ELISA Kit
Standar...