


Product Datasheet

Human Tenascin-R ELISA Kit (orb546565)

Description	Human Tenascin-R ELISA Kit
Reactivity	Human
Range	156 pg/ml - 10,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: Q32-F1358
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum and plasma (heparin, EDTA).
Uniprot ID	Q92752
Sensitivity	50 pg/ml
Dilution Range	Dilution ratio of 1:1, concentration in serum and plasma is around 1200 pg/ml.
Expiration Date	Please enquire.



Human
Tenascin-R
PicoKine
ELISA Kit
Stan...