


Product Datasheet

Human ANGPTL6 ELISA Kit (orb546574)

Description	Human ANGPTL6 ELISA Kit
Reactivity	Human
Range	31.2 pg/ml - 2,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: CHO; Immunogen sequence: Ala20-Leu470
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum and plasma (heparin or EDTA)
Uniprot ID	Q8NI99
Sensitivity	10 pg/ml
Dilution Range	Dilution ratio of 1:1, concentration in serum and plasma is around 300 pg/ml.
Expiration Date	Please enquire.



Human
ANGPTL6
PicoKine
ELISA Kit
Standar...