


Product Datasheet

Rat SCGB3A1 ELISA Kit (orb546575)

Description	Rat SCGB3A1 ELISA Kit
Reactivity	Rat
Range	15.6 pg/ml - 1,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: E.coli; Immunogen sequence: F22-G104
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum and plasma (heparin or EDTA)
Uniprot ID	E9PU36
Sensitivity	10 pg/ml
Dilution Range	Dilution ratio of 1:1, concentration in serum and plasma is around 60 pg/ml.
Expiration Date	Please enquire.



Rat
SCGB3A1
PicoKine
ELISA Kit
Standard ...