


Product Datasheet

Human PSG1 ELISA Kit (orb546578)

Description	Human PSG1 ELISA Kit
Reactivity	Human
Range	1.56 ng/ml - 100 ng/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: Q35-P419
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum, plasma (heparin or EDTA) and urine
Uniprot ID	P11464
Sensitivity	50 pg/ml
Dilution Range	Dilution ratio of 1:100, concentration in serum and plasma is around 1.25 ug/ml.
Expiration Date	Please enquire.



Human PSG1
PicoKine
ELISA Kit
Standard C...