


Product Datasheet

Mouse LI-Cadherin-17 CDH17 ELISA Kit (orb546579)

Description	Mouse LI-Cadherin-17 CDH17 ELISA Kit
Reactivity	Mouse
Range	46.9 pg/ml - 3,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: F26-S827
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum and plasma (heparin)
Uniprot ID	Q9R100
Sensitivity	10 pg/ml
Dilution Range	Dilution ratio of 1:400, concentration in serum and plasma is around 150 ng/ml.
Expiration Date	Please enquire.



Mouse
Cadherin-
17 PicoKine
ELISA Kit
Sta...