


Product Datasheet

Human Nephtrin ELISA Kit (orb546580)

Description	Human Nephtrin ELISA Kit
Reactivity	Human
Range	156 pg/ml - 10,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: Q23-S1055
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum and plasma (heparin, EDTA or citrate)
Uniprot ID	O60500
Sensitivity	Less than 50 pg/ml
Dilution Range	Dilution ratio of 1:1, concentration in serum and plasma is around 200 pg/ml.
Expiration Date	Please enquire.



Human
Nephtrin
PicoKine
ELISA Kit
Standar...