


Product Datasheet

Human LTBR ELISA Kit (orb546582)

Description	Human LTBR ELISA Kit
Reactivity	Human
Range	15.6 pg/ml - 1,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: CHO; Immunogen sequence: Q31-M227
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum and plasma (heparin, EDTA).
Uniprot ID	P36941
Sensitivity	10 pg/ml
Dilution Range	For serum, dilution ratio is 1:1 and concentration is around 300 pg/ml. For plasma, dilution ratio is 1:1 and concentration is around 60 pg/ml.
Expiration Date	Please enquire.



Human LTBR
PicoKine
ELISA Kit
Standard C...