


Product Datasheet

Mouse Matrilin-2 ELISA Kit (orb546601)

Description	Mouse Matrilin-2 ELISA Kit
Reactivity	Mouse
Range	93.7 pg/ml - 6,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: R24-R956
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum and plasma (heparin, EDTA).
Uniprot ID	008746
Sensitivity	50 pg/ml
Dilution Range	Dilution ratio of 1:5, concentration in serum and plasma is 300-3200 pg/ml.
Expiration Date	Please enquire.



Mouse
Matrilin-2
PicoKine
ELISA Kit
Stan...