


Product Datasheet

Human SMPD1 ELISA Kit (orb546605)

Description	Human SMPD1 ELISA Kit
Reactivity	Human
Range	156 pg/ml - 10,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: CHO; Immunogen sequence: L47-C631
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum and plasma (heparin).
Uniprot ID	P17405
Sensitivity	10 pg/ml
Dilution Range	Dilution ratio of 1:4, concentration in serum and plasma is around 6 ng/ml.
Expiration Date	Please enquire.



Human
SMPD1
PicoKine
ELISA Kit
Standard ...