


Product Datasheet

Mouse ERAP1 ELISA Kit (orb546608)

Description	Mouse ERAP1 ELISA Kit
Reactivity	Mouse
Range	312 pg/ml - 20,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: M1-L930
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum and plasma (heparin, EDTA).
Uniprot ID	Q9EQH2
Sensitivity	50 pg/ml
Dilution Range	Dilution ratio of 1:100, concentration in serum and plasma is 0.25-1.3 ug/ml.
Expiration Date	Please enquire.



Mouse
ERAP1
PicoKine
ELISA Kit
Standard ...