


Product Datasheet

Mouse CHIT1 ELISA Kit (orb546613)

Description	Mouse CHIT1 ELISA Kit
Reactivity	Mouse
Range	31.2 pg/ml - 2,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: A22-S464
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum and plasma (heparin, EDTA).
Uniprot ID	Q9D7Q1
Sensitivity	50 pg/ml
Expiration Date	Please enquire.



Mouse CHIT1
PicoKine
ELISA Kit
Standard ...