

## Product Datasheet

# Human TCN2/Transcobalamin-2 ELISA Kit (orb546618)

**Description** Human TCN2/Transcobalamin-2 ELISA Kit

**Reactivity** Human

**Range** 156 pg/ml - 10,000 pg/ml

**Tested Applications** ELISA

**Immunogen** Expression system for standard: NS0; Immunogen sequence: E19-W427

**Storage** Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)

**Note** For research use only


**Sample Types** cell culture supernatants, serum, plasma (heparin, EDTA) and urine.

**Uniprot ID** **P20062**

**Sensitivity** 50 pg/ml

**Dilution Range** Dilution ratio of 1:10, concentration in serum and plasma is 18-22.4 ng/ml.

**Expiration Date** Please enquire.



Human  
TCN2  
PicoKine  
ELISA Kit  
Standard  
C...