


Product Datasheet

Human CLEC4C ELISA Kit (orb546626)

Description	Human CLEC4C ELISA Kit
Reactivity	Human
Range	15.6 pg/ml - 1,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: Phe46-Ile213
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum and plasma (heparin, EDTA, citrate).
Uniprot ID	Q8WTT0
Sensitivity	10 pg/ml
Expiration Date	Please enquire.



Human
CLEC4C
PicoKine
ELISA Kit
Standard...