


Product Datasheet

Human CLEC10A ELISA Kit (orb546628)

Description	Human CLEC10A ELISA Kit
Reactivity	Human
Range	15.6 pg/ml - 1,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: Q61-H316
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum and plasma (heparin, EDTA, citrate).
Uniprot ID	Q8IUN9
Sensitivity	10 pg/ml
Dilution Range	Dilution ratio of 1:1, concentration in serum and plasma is 0-40 pg/ml.
Expiration Date	Please enquire.



Human
CLEC10A
PicoKine
ELISA Kit
Standar...