


Product Datasheet

Mouse SEMA4G ELISA Kit (orb546634)

Description	Mouse SEMA4G ELISA Kit
Reactivity	Mouse
Range	62.5 pg/ml - 4,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: V18-M673
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum and plasma (heparin, EDTA, citrate).
Uniprot ID	Q9WUH7
Sensitivity	10 pg/ml
Dilution Range	Dilution ratio of 1:1, concentration in serum and plasma is around 300 pg/ml.
Expiration Date	Please enquire.



Mouse
SEMA4G
PicoKine
ELISA Kit
Standard...