

Product Datasheet

Mouse SEMA6A/Semaphorin-6A ELISA Kit (orb546636)

Description Mouse SEMA6A/Semaphorin-6A ELISA Kit

Reactivity Mouse

Range 62.5 pg/ml - 4,000 pg/ml

Tested Applications ELISA

Immunogen Expression system for standard: NS0; Immunogen sequence: G19-T649

Storage Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)


Note For research use only

Sample Types cell culture supernatants, serum and plasma (heparin, EDTA, citrate).

Uniprot ID **O35464**

Sensitivity 50 pg/ml

Expiration Date Please enquire.



Mouse
SEMA6A
PicoKine
ELISA Kit
Standard...