


Product Datasheet

Human SEMA6D ELISA Kit (orb546640)

Description	Human SEMA6D ELISA Kit
Reactivity	Human
Range	31.2 pg/ml - 2,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: V21-N662
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum and plasma (heparin, EDTA, citrate).
Uniprot ID	Q8NFY4
Sensitivity	20 pg/ml
Dilution Range	Dilution ratio of 1:1, concentration in serum and plasma is less than the lowest standard, 31.2 pg/ml.
Expiration Date	Please enquire.



Human
SEMA6D
PicoKine
ELISA Kit
Standard...