


Product Datasheet

Mouse RTN4R ELISA Kit (orb546646)

Description	Mouse RTN4R ELISA Kit
Reactivity	Mouse
Range	156 pg/ml - 10,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: C27-S447
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum and plasma (heparin, EDTA, citrate).
Uniprot ID	Q99PI8
Sensitivity	50 pg/ml
Expiration Date	Please enquire.



Mouse
RTN4R
PicoKine
ELISA Kit
Standard ...