


Product Datasheet

Mouse BR-Cadherin-12 Cdh12 ELISA Kit (orb546648)

Description	Mouse BR-Cadherin-12 Cdh12 ELISA Kit
Reactivity	Mouse
Range	156 pg/ml - 10,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: L18-A716
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum and plasma (heparin, EDTA).
Uniprot ID	Q5RJH3
Sensitivity	50 pg/ml
Dilution Range	Dilution ratio of 1:50, concentration in serum and plasma is 60-130 ng/ml.
Expiration Date	Please enquire.



Mouse
Cadherin-
12 PicoKine
ELISA Kit
Sta...