


Product Datasheet

Human GFRA1 ELISA Kit (orb546649)

Description	Human GFRA1 ELISA Kit
Reactivity	Human
Range	46.9 pg/ml - 3,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: D25-S429
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum and plasma (heparin, EDTA, citrate).
Uniprot ID	P56159
Sensitivity	20 pg/ml
Expiration Date	Please enquire.



Human
GFRA1
PicoKine
ELISA Kit
Standard ...