


Product Datasheet

Mouse CD99 ELISA Kit (orb546669)

Description	Mouse CD99 ELISA Kit
Reactivity	Mouse
Range	125 pg/ml - 8,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: E.coli; Immunogen sequence: D29-G137
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum and plasma (heparin, EDTA, citrate).
Uniprot ID	Q8VCN6
Sensitivity	50 pg/ml
Expiration Date	Please enquire.



Mouse CD99
PicoKine
ELISA Kit
Standard C...