


Product Datasheet

Human M-Cadherin-15 CDH15 ELISA Kit (orb546679)

Description	Human M-Cadherin-15 CDH15 ELISA Kit
Reactivity	Human
Range	156 pg/ml - 10,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: L61-A606
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum and plasma (heparin, EDTA, citrate).
Uniprot ID	P55291
Sensitivity	50 pg/ml
Expiration Date	Please enquire.



Human
Cadherin-
15 PicoKine
ELISA Kit
Sta...