

Product Datasheet

Human 5-HT2B Receptor protein (orb54671)

Description	Human 5-HT2B Receptor protein
Endotoxins	Not test.
Preservatives	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.
Form/Appearance	Liquid or Lyophilized powder
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Tag	N-terminal 6xHis-SUMO-tagged
Note	For research use only
Application notes	N-terminal 6xHis-SUMO-tagged: N-terminal 6xHis-SUMO-tagged1-352AA: 1-481AA Full Length : Full Length
Protein Sequence	MALSYRVSELSQSTPEHILQSTFVHVHSSNWSGLQTESIPEEMKQIVVEEQGNKLHW AALLILMVIPTIGGNTLVILAVSLEKKLQYATNYFLMSLAVADLLVGLFVMPIALTI MFEAMWPLPLVLCPAWLFLDVLFTASIMHLCAISVDYIAIKKPIQANQYNSRAT AFIKITVWVLISIGIAIPVPIKGIETDNDPNNITCVLTKERFGDFMLFGSLAAFFTPL AIMIVTYFLTIHALQKKAYLVKNKPPQRLTWLTVSTVFQRDETPCSSPEKVAMLDG SRKDKALPNSGDETLMRRTSTIGKKSQVTSNEQRASKVLGIVFFLFLMWCPFFI TNITLVLCDSNQTTLQMLLEIFVWIGYVSSGVNPLVYTLFNKTFRDAFGRYITCN YRATKSVKTLRKSSKIYFRNPMAENSKFFKKHGIRNGINPAMYQSPMRLRSSTIQ SSSIILLDTLLLTENEGDKTEEQVSYV
Purity	Greater than 90% as determined by SDS-PAGE.
Source	in vitro E.coli expression system
MW	70.3 kDa
Uniprot ID	P41595
Expiration Date	6 months from date of receipt.