


Product Datasheet

Mouse sE-Selectin ELISA Kit (orb546909)

Description	Mouse sE-Selectin ELISA Kit
Reactivity	Mouse
Range	156 pg/ml - 10,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: W22-P557
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum and plasma (heparin, citrate).
MW	66750 MW
Uniprot ID	Q00690
Sensitivity	10 pg/ml
Expiration Date	Please enquire.



Mouse sE-
Selectin
PicoKine
ELISA Kit
sta...