


Product Datasheet

Mouse sL-Selectin ELISA Kit (orb546911)

Description	Mouse sL-Selectin ELISA Kit
Reactivity	Mouse
Range	62.5 pg/ml - 4,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: W39-N332
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
MW	38463 MW
Uniprot ID	Q3TCF3
Sensitivity	5 pg/ml
Expiration Date	Please enquire.



Mouse L-
Selectin
PicoKine
ELISA Kit
stan...