

Product Datasheet

Mouse P-Selectin / CD62P ELISA Kit (orb546913)

Description	Mouse P-Selectin / CD62P ELISA Kit
Reactivity	Mouse
Range	312 pg/ml - 20,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: W42-A709
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum and plasma (citrate)
MW	83099 MW
Uniprot ID	Q01102
Sensitivity	5 pg/ml
Expiration Date	Please enquire.

