

Product Datasheet

Mouse VEGF ELISA Kit (orb546916)

Description	Mouse VEGF ELISA Kit
Reactivity	Mouse
Range	15.6 pg/ml - 1,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: A27-R190
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum and plasma (heparin, EDTA, citrate)
MW	25283 MW
Uniprot ID	Q00731
Sensitivity	2 pg/ml
Expiration Date	Please enquire.

