


Product Datasheet

Mouse CXCL5/LIX/ENA-78/GCP 2 ELISA Kit (orb546938)

Description	Mouse CXCL5/LIX/ENA-78/GCP 2 ELISA Kit
Reactivity	Mouse
Range	15.6 pg/ml - 1,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: E.coli; Immunogen sequence: V45-A118
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum and plasma (heparin, EDTA, citrate)
MW	14235 MW
Uniprot ID	P50228
Sensitivity	10 pg/ml
Expiration Date	Please enquire.



Mouse
CXCL5
PicoKine
ELISA Kit
standard ...