


Product Datasheet

Human Granulysin ELISA Kit (orb546949)

Description	Human Granulysin ELISA Kit
Reactivity	Human
Range	31.2 pg/ml - 2,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: M1-L145
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
MW	16374 MW
Uniprot ID	P22749
Sensitivity	10 pg/ml
Expiration Date	Please enquire.



Human
Granulysin
PicoKine
ELISA Kit
stan...