


Product Datasheet

Human CD21/CR2 ELISA Kit (orb546951)

Description	Human CD21/CR2 ELISA Kit
Reactivity	Human
Range	31.2 pg/ml - 2,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: CHO; Immunogen sequence: Ile21-71
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
MW	112916 MW
Uniprot ID	P20023
Sensitivity	20 pg/ml
Expiration Date	Please enquire.



Human
CD21
PicoKine
ELISA Kit
standard c...