


Product Datasheet

Human TGFBR3/Tgf Beta Riii ELISA Kit (orb546956)

Description	Human TGFBR3/Tgf Beta Riii ELISA Kit
Reactivity	Human
Range	156 pg/ml - 10,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: G21-D781
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
MW	93499 MW
Uniprot ID	Q03167
Sensitivity	10 pg/ml
Expiration Date	Please enquire.



Human
TGFBR3
PicoKine
ELISA Kit
standard...