

Product Datasheet

GFP Rabbit Monoclonal Antibody, HRP Conjugated (orb546965)

Description	GFP Rabbit Monoclonal Antibody, HRP Conjugated
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A synthesized peptide
Form/Appearance	Liquid
Concentration	Actual concentration vary by lot. Use suggested dilution ratio to decide dilution procedure.
Storage	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Monoclonal
MW	25993 MW
Uniprot ID	Q963I9
Dilution Range	WB 1:3000-1:10000
Expiration Date	12 months from date of receipt.

