



Product Datasheet

p16 INK Monoclonal Antibody (orb546982)


Description	p16 INK Monoclonal Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, ICC, IF, IHC, IP, WB
Immunogen	A synthesized peptide derived from human p16 INK
Form/Appearance	Liquid
Concentration	Actual concentration vary by lot. Use suggested dilution ratio to decide dilution procedure.
Storage	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Monoclonal
Uniprot ID	P42771
Dilution Range	WB 1:500-1:2000IHC 1:50-1:200ICC/IF 1:50-1:200IP 1:30FC 1:100
Expiration Date	12 months from date of receipt.



WB analysis of CDKN2A using anti-CDKN2A ...



IF analysis of immunocytochemical section...



IHC analysis of CDKN2A using anti-CDKN2A...