

Product Datasheet

pan Cadherin 1 CDH1 Rabbit Monoclonal Antibody (orb547021)

Description pan Cadherin 1 CDH1 Rabbit Monoclonal Antibody

Species/Host Rabbit

Reactivity Human, Mouse, Rat

Conjugation Unconjugated

Tested Applications FC, ICC, IF, IHC, IP, WB

Immunogen A synthesized peptide derived from human pan Cadherin

Form/Appearance Liquid

Concentration Actual concentration vary by lot. Use suggested dilution ratio to decide dilution procedure.

Storage Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.

Note For research use only

Isotype Rabbit IgG


Clonality Monoclonal

MW 97456 MW


Uniprot ID **P39038**, **P55283**, **P19022**, **P22223**, **P33151**, **P12830**

Dilution Range WB 1:1000-1:5000IHC 1:50-1:200ICC/IF 1:50-1:200IP 1:50FC 1:50

Expiration Date 12 months from date of receipt.



WB analysis of pan Cadherin expression i...



Immunohistochemical analysis of paraffin...