

Product Datasheet

HMGB1/Hmg 1 Rabbit Monoclonal Antibody (orb547024)

Description HMGB1/Hmg 1 Rabbit Monoclonal Antibody

Species/Host Rabbit

Reactivity Human, Mouse, Rat

Conjugation Unconjugated

Tested Applications FC, ICC, IF, IHC, WB

Immunogen A synthesized peptide derived from human HMGB1

Form/Appearance Liquid

Concentration Actual concentration vary by lot. Use suggested dilution ratio to decide dilution procedure.

Storage Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.

Note For research use only

Isotype Rabbit IgG


Clonality Monoclonal

MW 24894 MW


Uniprot ID **P09429**

Dilution Range WB 1:1000-1:2000IHC 1:50-1:100ICC/IF 1:50-1:100FC 1:30

Expiration Date 12 months from date of receipt.



WB analysis of
HMGB1 expression in
(1)He...



Immunohistochemical
analysis of paraffin...