


Product Datasheet

PCSK9/Proprotein Convertase 9 Rabbit Monoclonal Antibody (orb547042)

Description	PCSK9/Proprotein Convertase 9 Rabbit Monoclonal Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IP, WB
Immunogen	A synthesized peptide derived from human PCSK9
Form/Appearance	Liquid
Concentration	Actual concentration vary by lot. Use suggested dilution ratio to decide dilution procedure.
Storage	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Monoclonal
MW	74286 MW
Uniprot ID	Q8NBP7
Dilution Range	WB 1:500-1:2000IP 1:20
Expiration Date	12 months from date of receipt.



WB analysis of PCSK9 expression in HepG2...